

# A new concept for isotope ratio monitoring LC/MS

## A Wide Range of Applications

**Andreas Hilkert, Dieter Juchelka, Michael Krummen**

Thermo Electron (Bremen) GmbH

# Overview


- Introduction in *Isotope Ratio Monitoring*-LC/MS (*irm*-LC/MS)
  - *Technology*
  - *Operating Modes*
- New Applications by *irm*-LC/MS
  - *Authenticity Control*
    - Detection of adulteration by sugars in honey.
  - *Determination of Origin*
    - Differentiation of analgesic drugs.
  - *Molecular Biology*
    - Carbon isotopic characterization of rRNA.
  - *Biogeochemistry*
    - Plant metabolism study of organic acids.
  - *Forensic Chemistry*
    - Analysis of aspartic acid in cadaver blood samples.

# Why irm – LC/MS ?

- $\delta^{13}\text{C}$  analysis of individual compounds with:
  - *High molecular weight*
  - *High polarity*
  - *Thermal instability*
  - *Low vapour pressure*
  
  - *Less sample preparation*
  - *No derivatization*
  - *No isotope dilution*
  - *Less risk of fractionation*

# Comparison between irm GC/MS and irm LC/MS

- irm – GC/MS




- Helium as carrier for
  - Separation of compounds
  - Transfer to the IRMS
- Helium has
  - No impact on combustion
  - No effects in the IRMS

**Dry combustion (oxidation)  
in the He phase**

# Comparison between irm GC/MS and irm LC/MS

- irm-LC/MS (first strategy)




- Solvents as carrier for
  - Separation of compounds
- Solvents are
  - Oxidized
  - Hazardous to IRMS
- No solvents to reactor or IRMS


## First approaches ('91, '93)

- “Moving Wire” – Drying system (difficult to use)
- “Particle Beam” Separation (low sensitivity, fractionation)

# A New Strategy




# Scheme of the LC IsoLink Interface



# HPLC Resolution

➔ HPLC separation of carbohydrates




➔ HPLC resolution is maintained

- Parameters:
  - HPLC flow:
    - 300  $\mu\text{l}/\text{min}$
  - Oxidation reagent:
    - 60  $\mu\text{l}/\text{min}$   
( $\text{NH}_4$ )<sub>2</sub>S<sub>2</sub>O<sub>8</sub>, 100g/l
  - Column:
    - 700 CH  
Carbohydrate  
Column, 90 °C
  - Reactor:
    - 99.9 °C
  - CO<sub>2</sub> Exchanger:
    - 1 ml/min He flow



# Authenticity Control of Honey

Investigation of the adulteration of honey analyzing glucose, sucrose, fructose




Honey	Glucose $\delta^{13}\text{C}\text{‰}$	Fructose $\delta^{13}\text{C}\text{‰}$	Area Fru/Glu	
A	-27.9	-27.8	1.13	pure
B	-25.1	<b>-26.4</b>	2.17	<b>adulterated</b>
C	-26.5	-26.5	1.35	pure
D	-26.1	-26.0	<b>4.53</b>	<b>adulterated</b>
E	-11.2	<b>-13.9</b>	0.65	<b>adulterated</b>

- Absolute  $\delta^{13}\text{C}$  value
- $\delta^{13}\text{C}$  difference, Glu – Fru
- Ratio of area, Fru / Glu

# Source Differentiation of Drugs


➔ Determination of different analgesic compounds.  
 $\delta^{13}\text{C}$  of Paracetamol (Acetamidophenol) and Aspirin<sup>®</sup> (Acetylsalicylic Acid; ASA).



Tablet Type	$\delta^{13}\text{C}$ (‰)		
	Paracetamol	ASA	$\mu$ -EA
A Country1	-	<b>-34.2</b>	<b>-33.4</b>
A Country2	-	<b>-34.2</b>	<b>-33.5</b>
B	-	-29.1	-27.6
C	-	<b>-27.2</b>	<b>-26.6</b>
D <sub>1</sub>	-	<b>-26.8</b>	<b>-26.7</b>
E	-	-27.7	-27.6
F	-32.3	-32.6	-31.7
D <sub>2</sub>	-28.7	<b>-33.8</b>	-31.3
G	-29.2	-32.7	-29.7

- Tablet type A has the same origin
- 4 sources of ASA
- Producer D use different ASA sources


# $\mu$ -EA – Direct Injector



- Fast analysis of all water soluble compounds

HPLC Analysis

$\mu$ -EA




Analysis of a tablet followed by direct loop injection ( $\mu$ -EA). Loop size of the HPLC injector was 5  $\mu$ L, the loop size of the  $\mu$ -EA injector was 10  $\mu$ L, which results in two-fold response of the  $\mu$ -EA peak

# μ-EA – Reproducibility

## ➔ Bulk Injection of Amphetamine

m/z 44 (mV)




e.g., 5 x bulk injections of 218 ng amphetamine


Sample	Amount		$\delta^{13}\text{C}$ (‰)	S.D. (‰)	n
	Amphetamine (ng)	Carbon (ng)			
#1	218	174	-28.70	0.04	5
	452	362	-28.76	0.04	5
	698	558	-28.83	0.03	5
	944	756	-28.87	0.03	5
	1161	929	-28.76	0.03	5
	1359	1088	-28.67	0.02	5
	2280	1824	-28.70	0.02	5
<b>Mean</b>			<b>-28.76</b>	<b>0.03</b>	
#2	732	586	<b>-32.23</b>	<b>0.04</b>	5
#3	438	350	<b>-31.58</b>	<b>0.04</b>	5
#4	781	625	<b>-27.89</b>	<b>0.03</b>	5

➔ Reliable reproducibility of the  $\delta^{13}\text{C}$  values

# irm-LC/MS of NaOH-hydrolyzed *E. coli* RNA




Linkage of microbiological species identity with carbon source utilization by carbon isotopic characterization of rRNA.



A: Adenosine Nucleotides


- Cell growth and Isolation: RNA was extracted from overnight cultures of *E. coli* grown on M9 minimal salts with glucose as sole carbon source Medium.
- Hydrolyse: 0.2 N NaOH for 15 min at 50 °C

# Carbon Isotopic Composition of Individual Peaks from NaOH-hydrolyzed E. coli RNA




⇒ The carbon isotopic composition of RNA closely reflects that of the growth substrate.

# $\delta^{13}\text{C}$ Analysis of Fruit Juice Organic Acids



- **HPLC Column:**  
Allure Organic Acids,  
300 mm x 4.6 mm, 5  $\mu\text{m}$
- **Flow:**  
500  $\mu\text{l}/\text{min}$
- **Mobile Phase:**  
100 mM  $\text{KH}_2\text{PO}_4$ , pH 3


# $\delta^{13}\text{C}$ Analysis of Organic Acids






# *irm*-LC/MS: $\delta^{13}\text{C}$ Analysis of an Extract of *Geranium pratense* (May 2004)

## Plant metabolism study of organic acids



# Analysis of Volatile Fatty Acids by irm-LC/MS



1. **Formate**
2. **Lactate**
3. **Acetate**
4. **Propionate**

Pump 1: 200 mg/l  $(\text{NH}_4)_2\text{S}_2\text{O}_8$ ,  
50  $\mu\text{l}/\text{min}$

Pump 2: 2M  $\text{H}_3\text{PO}_4$ ,  
no  $\text{AgNO}_3$ -catalyst,  
50  $\mu\text{l}/\text{min}$


Reactor: 99.9°C

Data: Prof. Hirrichs Univ. Bremen


# Linearity of the IRMS - Signal

Amounts injected : 8.8 nmol – 0.28 nmol


Formate




Lactate



Acetate




Propionate



Data: Prof. Hirrichs Univ. Bremen

# $\delta^{13}\text{C}$ Analysis of Aspartic Acid in Cadaver Blood



<b>Aspartic Acid</b>		
Sample	$\delta^{13}\text{C}$ (‰)	Std. Dev.
A	-32.12	0.08
B	-32.66	0.11
C	-16.18	0.03
D	-18.76	0.14
E	-19.98	0.06
F	-16.62	0.12
Standard		
G	-24.79	0.05
H	-26.19	0.11

- Mobile Phase: 10mM  $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$ , pH 4.7, Flow: 300  $\mu\text{l}/\text{min}$
- HPLC Column: Nova-Pack<sup>®</sup> C18, 60 Å (4  $\mu\text{m}$ , 3.9 mm x 300 mm).

# Summary

- *irm*-LC/MS opens a wide range of interesting applications.
- Macromolecules, non-volatile components and components which tend to decompose are directly accessible for precise isotopic analysis.

# Acknowledgements

- Barbara J. MacGregor and Adam Friedman, Department of Marine Sciences, University of North Carolina, Chapel Hill, N.C.
- Elena Hettmann and Gerd Gleixner, Max-Planck-Institut für Biogeochemie, Jena
- Jean-Philippe Godin and Jörg Hau, Nestlé Research Center, Lausanne
- Kurt-Peter Raezke, APPLICA GmbH, Bremen